

Table I. *Semen collection.* Two months after surgery, the semen was proved to be of good quality both in terms of concentration and vitality.

Subject	1° collection	2° collection
Reaction time (min.)	10 minutes	15 minutes
Volume before filtration (ml)	100	80
Volume after filtration (ml)	85	70
Color	Milky	Milky
Smell	Heated metal	Heated metal
pH	7.3	7.3
Motility (%)	85	90
Vitality (%)	90	95
Concentration (X 10 ⁶ /ml)	210	145
Total number of spermatozoa for ejaculate (X 10 ⁹)	17	10.150
Motility at environmental temperature (h)	2	2.30
Primary morphological abnormalities (%)	10	7
Secondary morphological abnormalities (%)	12	9